

05462	cagar	1KI -6: 20	pure
05462	cagar	1KI -6: 21	pure
05462	cagar	1KI -7: 49	pure
05462	cagar	1KI -7: 50	pure
05462	cagar	1KI -10: 21	pure
05462	cagar	2CH-4: 20	pure
05462	cagar	2CH-4: 22	pure
05462	cagar	2CH-9: 20	pure
05462	cagar	1KI -11: 27	repa red
05462	cagar	GEN-7: 16	shut
05462	cagar	GEN-19: 6	shut
05462	cagar	GEN-19: 10	shut
05462	cagar	EX0-14: 3	shut
05462	cagar	LEV-13: 4	shut
05462	cagar	LEV-13: 5	shut
05462	cagar	LEV-13: 11	shut
05462	cagar	LEV-13: 21	shut
05462	cagar	LEV-13: 26	shut
05462	cagar	LEV-13: 31	shut
05462	cagar	LEV-13: 33	shut
05462	cagar	LEV-13: 50	shut
05462	cagar	LEV-13: 54	shut
05462	cagar	LEV-14: 38	shut
05462	cagar	LEV-14: 46	shut
05462	cagar	NUM-12: 14	shut
05462	cagar	NUM-12: 15	shut
05462	cagar	DEU-32: 30	shut
05462	cagar	JOS-2: 7	shut
05462	cagar	JOS-6: 1	shut
05462	cagar	JUD-3: 23	shut
05462	cagar	JUD-9: 51	shut
05462	cagar	1SA-1: 5	shut
05462	cagar	1SA-1: 6	shut

05462	cagar	1SA-23: 7	shut
05462	cagar	2KI -4: 4	shut
05462	cagar	2KI -4: 5	shut
05462	cagar	2KI -4: 21	shut
05462	cagar	2KI -4: 33	shut
05462	cagar	2KI -6: 32	shut
05462	cagar	2CH-28: 24	shut
05462	cagar	2CH-29: 7	shut
05462	cagar	NEH-6: 10	shut
05462	cagar	NEH-13: 19	shut
05462	cagar	JOB-3: 10	shut
05462	cagar	JOB-11: 10	shut
05462	cagar	JOB-41: 15	shut
05462	cagar	PSA-31: 8	shut
05462	cagar	ECC-12: 4	shut
05462	cagar	I SA-22: 22	shut
05462	cagar	I SA-22: 22	shut
05462	cagar	I SA-24: 10	shut
05462	cagar	I SA-24: 22	shut
05462	cagar	I SA-26: 20	shut
05462	cagar	I SA-45: 1	shut
05462	cagar	I SA-60: 11	shut
05462	cagar	JER-13: 19	shut
05462	cagar	EZE-3: 24	shut
05462	cagar	EZE-44: 1	shut
05462	cagar	EZE-44: 2	shut
05462	cagar	EZE-44: 2	shut
05462	cagar	EZE-46: 1	shut
05462	cagar	EZE-46: 2	shut
05462	cagar	EZE-46: 12	shut
05462	cagar	MAL-1: 10	shut
05462	cagar	JOB-12: 14	shutteth
05462	cagar	JOS-2: 5	shutting

05462 cagar PSA-35: 3 stop  
05462 cagar PSA-78: 48 gave  
05462 cagar PSA-78: 50 gave  
05462 cagar PSA-78: 62 gave  
05462 cagar LAM-2: 7 gi ven  
05462 cagar PSA-17: 10 i ncl osed  
05462 cagar PSA-78: 50 over  
05462 cagar PSA-78: 62 over  
05462 cagar GEN-2: 21 cl osed  
05462 cagar JUD-3: 22 cl osed  
05462 cagar DEU-23: 15 del i ver  
05462 cagar JOS-20: 5 del i ver  
05462 cagar 1SA-17: 46 del i ver  
05462 cagar 1SA-23: 11 del i ver  
05462 cagar 1SA-23: 12 del i ver  
05462 cagar 1SA-23: 12 del i ver  
05462 cagar 1SA-23: 20 del i ver  
05462 cagar 1SA-30: 15 del i ver  
05462 cagar AMO-1: 6 del i ver  
05462 cagar AMO-6: 8 del i ver  
05462 cagar 1SA-24: 18 del i vered  
05462 cagar 1SA-26: 8 del i vered  
05462 cagar 2SA-18: 28 del i vered  
05462 cagar JOB-16: 11 del i vered  
05462 cagar AMO-1: 9 del i vered  
05462 cagar OBA-1: 14 del i vered

~~~~~

~~~~~